

Observations at a PFAS Contaminated Site: Variability and Precursor Occurrence

Carolyn Acheson, John Carriger, Sandra Halstead, Ronald Herrmann, Danielle Kleinmaier, Marc Mills, and Lawrence Zintek (US EPA)

and

Solidea M.C. Bonina, Andrea Burkes, and Raghuraman Venkatapathy (Pegasus Technical Services, Inc.)

Office of Research and Development National Risk Management Research Laboratory Land and Materials Management Division



Disclaimer

The views expressed in this presentation are those of the authors and do not necessarily reflect the views or policies of the U.S. Environmental Protection Agency.



Joint Base Elmendorf Richardson

- Located in Anchorage, AK
- Active from early 1940's
- Listed on NPL
- 14 OUs+ combined
 - Size 86,000 acres
 - Geology coastal, glacial, forested, mountainous
 - Contaminants solvents, landfills, munitions
- Air Force (AF) lead in 2010
- Anchorage drinks surface water





2016 JBER Sampling Effort

- 26 Areas of Concern investigated
 - Flightline, hangars, fire stations, fire training pits, crash sites
 - Media included ground water and soils
 - Later ground water seeps and surface water collected

EPA

- EPA collected samples immediately after AF collected samples
- 17 ground water locations
- 4 auxiliary locations not sampled by AF
- 6 seeps sampled later





EPA objectives for JBER sampling

- Evaluate EPA analytical method on new matrices
 - Accuracy
 - Precision
- Evaluate sample variability in replicate samples
 - 3 field replicates for many locations
- Analyze samples for a larger suite of PFAS including precursors and transformation products
 - 12 PFCAs: C4 C14
 - 7 PFSAs: C4 C10
 - 12 precursors





EPA used ASTM Method D7979

- Environmental Waters (non-potable)
- Direct Injection, analysis by LC/MS/MS
- Single laboratory validated
- Target Analytes:
 - 12 PFCAs C4 C14
 - 7 PFSAs C4 C10
 - 12 precursors
- Isotopically labeled surrogates:
 - 7 PFCAs, 2 PFSAs, 3 precursors
 - Used to monitor analytical method performance/quality
 - Not used to "correct" the data
- Uses confirmation ion ratios to identify compounds and minimize matrix issues
- ASTM D7979 updated since this study was conducted





ASTM D7979 method



Based on schematic by William Lipps, Shimadzu



Analytical Method Quality Controls

- Analyte Identification
 - Each batch: Initial calibration, Calibration check, and Second source check
 - Each analyte: Retention time, Primary and Confirmation ion masses, and Ion ratio
- Accuracy 2 of each/batch unless specified
 - Surrogate spiking All samples and blanks
 - Used to assess method performance
 - Not used to alter reported concentrations
 - Matrix spike samples MS and MS duplicates
 - Spiked blanks
 - Method reporting limit checks
- Precision 2 of each/batch
 - Duplicate samples
 - Matrix spike duplicates
 - Spiked blanks

7

Laboratory Contamination – method blanks – 2/batch





Deviations from ASTM D7979 in JBER sampling and analysis

- 50 mL sample collected required sample processing changes (purple arrows)
- Sample coolers arrived above 6 C in some cases
- Sample holding time exceeded





Most commonly observed **PFAS in JBER samples**



- N= 77 samples
- **Reporting level**
 - 10 ng/L most
 - 15 ng/L PFOS
 - 50 ng/L PFPeA
 - 30 ng/L PFBS
 - Screening level
 - PFOA 70 ng/L
 - PFOS 70 ng/L •
 - PFBS 380 µg/L

9



Concentration ranges observed





Method Performance Accuracy - Matrix Spike Data





Method Performance Accuracy - Surrogate Recovery



SERA United States Environmental Protection Agency Matrix Spike and Lab Duplicates



N dup = 9 to 12 N MS Dup = 16

EPA
United States
Environmental ProtectionMethod Performance - PrecisionMatrix Spike and Lab Duplicates



N dup = 6 to 10 N MS Dup = 16



Additional Analytical Results

- Field and Equipment Blanks no analytes observed above reporting levels
- Performance Evaluation Sample
 - Double Blind sample
 - Results within duplicate data acceptance criteria

Chemical	Measured Conc (ng/L)	Spiked Conc (ng/L)	RPD (%)
PFOA	115	100	15
Surrogate recovery	99.6 %		
PFOS	210	200	5.0
Surrogate recovery	101 %		



Variability in Field Replicates





Variability in Field Replicates compared to Lab Replicates



RPD for field dup based on median

Field dup data not meeting DQOs are not shown

*1 outlier at 125% not shown

Field dups N= 37 to 42



Precursor observations

- Observed in 70 % of sample locations
- Analytes observed
 - 6:2 FTS 49 samples
 - 4:2 and 8:2 FTS 9 samples
 - 6:2 FTUCA and FOSA 1 sample
 - Not observed 6:2 FTCA, 8:2 FTCA, 8:2 FTUCA, 7:3 FTCA, 10:2 FTCA, N-EtFOSAA, N-MeFOSAA



FTS Precursors compared to Total PFAS – Molar basis



Total PFAS (M)



JBER Conclusions

- Analytical method used by EPA
 - accurately and precisely measured concentrations for most analytes
 - 6:2 FTS analysis further method development useful
- Sample variability
 - Variability in field replicates similar to lab precision data
 - <u>+</u> 20% for many locations at this site
- Precursors
 - 6:2, 8:2 and 4:2 FTS most commonly observed precursors
 - Observed in 70 % of sample locations
 - Precursors relative to total measured PFAS molar basis
 - Common 10 %
 - Could be as high as 40 %





- Use exploratory data interpretation techniques such as:
 - Principal Component Analysis (PCA)
 - Hierarchical Cluster Analysis (HCA)
 - Bayesian Networks and machine learning

to evaluate the data such as identifying similarities in sample locations and PFAS patterns



Sample location clusters

- All PFAS and surrogate data analyzed using PCA and HCA
- Seven clusters were identified:

Cluster	Total PFAS Conc (M)	PFOS/PFOA ratio
Fire stn 7	1.4 x 10 ⁻⁷	6
Hangar 18	4.7 to 5.1 x 10 ⁻⁸	50
GW seep OU5SP-11 and seep WCSW-2	3.3 to 4.9 x 10 ⁻⁸	2
FSFS and Fire stn 1	1.8 to 2.1 x 10 ⁻⁸	2
GW seep OU5SP1 and 2, and pump stn OU5CP	5.3 to 7.1x 10 ⁻⁹	1
CHD3, Hangar 10, Hangar 8, and Hangar 6	3.1 to 9.0 x 10 ⁻⁹	1 to 4
All other sample locations	< 3.3 x 10 ⁻⁹	







Hangar 18 0

Fire stn 7

Seep OU5SP-11 Seep WCW-02

FSFSOFire stn 1O

CHD3

Hangar 10

Hangar 8

Hangar 6

Seep OU5SP-01 Seep OU5SP-02 Seep OU5CP-02



Bayesian Networks

Which analytes are correlated to detections of precursors?

Using unsupervised machine learning and Pearson's correlation coefficient

62 ETHIN

8:2 FTS

82 FTS

0.6487

0.65

0.3948

PENA



PFDA



Acknowledgements

- Reg 10: Kira Lynch, Karin Fedderson
- Reg 5: Robert Snyder, Sylvia Griffin, George Schupp, and Dennis Wesolowski (retired)
- ORD: John McKernan, Steve Musson, Michael Moeykens, Steven Jones, Jim Voit, Ed Barth, and Michelle Simon
- Senior Environmental Employee: Marilyn Dapper
- Ch2m: Jennifer Ulrich, Michael Landen, Alison Wieland, Kristen Stevens, and Morgan Bruno
- Anchorage Water and Wastewater Utility: Rob Gustafson
- Pegasus Technical Services, Inc: Robert Grosser

